



**IASLC 2025 World  
Conference on Lung Cancer**

**SEPTEMBER 6-9, 2025 | BARCELONA, SPAIN**

**wclc.iaslc.org**       **#WCLC25**

# Ultrasensitive MRD improves detection in resectable stage I NSCLC

Zachary Coyne<sup>1</sup>, David M Kurtz<sup>2</sup>, Sameena Khan<sup>1</sup>, Jamie Feng<sup>1</sup>, Thomas K Waddell<sup>3</sup>, Kazuhiro Yasufuku<sup>3</sup>, Najib Saffiedine<sup>4</sup>, Michael Ko<sup>5</sup>, Laura Donahoe<sup>3</sup>, Andrew Pierre<sup>3</sup>, Shaf Keshavjee<sup>3</sup>, Mary R Rabey<sup>1</sup>, Lisa W Le<sup>6</sup>, Alex Salvarrey<sup>3</sup>, Morgan Schroeder<sup>2</sup>, Giselle Zamora-Chavez<sup>2</sup>, Sierra L Stowell<sup>2</sup>, Hayley Warsinske<sup>2</sup>, Ash Alizadeh<sup>7</sup>, Maximilian Diehn<sup>2,7</sup>, Natasha B Leigh<sup>1</sup>

<sup>1</sup>Division of Medical Oncology, Princess Margaret Cancer Centre/University Health Network, Toronto, ON, Canada

<sup>2</sup>Foresight Diagnostics, Aurora, CO, USA

<sup>3</sup>Division of Thoracic Surgery, University Health Network, Toronto, ON, Canada

<sup>4</sup>Division of Thoracic Surgery, Michael Garron Hospital, Toronto, ON, Canada

<sup>5</sup>Division of Thoracic Surgery, Unity Health Centre, Toronto, ON, Canada


<sup>6</sup>Department of Biostatistics, University Health Network, Toronto, ON, Canada

<sup>7</sup>Department of Radiation Oncology, Stanford University, Stanford, CA, USA

CONQUERING LUNG AND OTHER THORACIC CANCERS WORLDWIDE IN THE 21ST CENTURY

# Suboptimal detection rates in stage I NSCLC using first-generation MRD assay

- Plasma ctDNA is a valuable tool for recurrence risk stratification in resected NSCLC
- First generation tumour-informed ctDNA MRD assays have suboptimal clinical sensitivity, particularly in stage I lung adenocarcinoma (LUAD) patients
- Our experience with a first generation MRD assay in clinical stage I NSCLC demonstrated a pre-op detection rate of 22.7% and a post-op detection rate of 2%<sup>1</sup>
- These findings highlight the need for **ultrasensitive ctDNA-MRD assays** to improve early detection and guide treatment decisions in stage I NSCLC

**Objectives:**

- Detection rate
- Sensitivity
- Association with RFS/OS

**NCT05254782**<sup>1</sup>Khan et al. ASCO 2024




# Analysis of banked samples from resected clinical stage I NSCLC with an ultrasensitive assay built on PhasED-Seq

## ctDNA-MRD Test Workflow

### IDENTIFICATION OF PATIENT-SPECIFIC VARIANTS




### MRD DETECTION



Analytical LOD95\* = 0.3 parts per million ( $3 \times 10^{-7}$ )

- This ctDNA-MRD test includes up to 5,000 patient-specific loci:
  - Phased variants (e.g., PhasED-Seq)
  - Additional alterations with low error profiles (e.g., selected SNVs) to maximize sensitivity




\*The LOD95 was determined in an analytical study using contrived cfDNA-into-cfDNA samples in a limited dilution series and modeled by Probit regression

Cabel et al. *ESMO* 2024

# Cohort demographics & methods


Characteristic		N=30 (%)
Clinical Stage	I	30 (100%)
Pathological Stage	IA	9 (30)
	IB	7 (23.3)
	II	8 (26.7)
	III	5 (16.7)
	IV	1 (3.3)
Histology	Adenocarcinoma	25 (83.3)
	Squamous	3 (10)
	Other	2 (6.7)
Smoking status	Never	12 (40)
	Former	15 (50)
	Current	3 (10)
Sex	Female	23 (76.7)
Median age (range)		67 (38-85)
Adjuvant therapy		17 (57.0)



- Panels designed starting from WGS and selecting up to 5000 targets, focused on phased variants
  - 100% success in panel design (30/30)
  - Median 4055 variants tracked
- cfDNA input for library preparation: median 18.5ng (8.5 - 48.0)
  - Plasma volume: median 4.9mL (3.8 – 6.0)

# Detection of post-surgical ctDNA MRD in NSCLC patients with relapse using ultrasensitive assay


- This cohort (n=30) was derived from the ctDNA-Lung-DETECT study which used a 1<sup>st</sup> generation MRD assay in patients with Stage I NSCLC reported at ASCO 2024<sup>1</sup>
- In this selected subgroup (n=30), the PhasED-Seq-based assay demonstrated ctDNA detection in Stage I NSCLC\*:
  - Pre-operative detection 68% (n=28)
  - Post-operative detection 38% (n=29)
  - Post-operative MRD detected in 55% of relapsed cases (n=11)
- Detectable MRD post-operatively was associated with adverse RFS
  - Post-operative landmark: HR = 3.14, P = 0.0425
  - 1-year post-op landmark: HR = 8.20, P = 0.0001



\* This is not a direct patient level comparison.

<sup>1</sup>Khan et al. ASCO 2024


# Post-operative ctDNA MRD detection in cases with relapse



- 11 cases with recurrence assessed
- In patients with MRD detected at the post-operative landmark, the median lead-time to recurrence was 10 months
- 11/12 (92%) of samples collected within 12 months preceding relapse had MRD detected

# Concluding remarks

- Ultrasensitive ctDNA-MRD detection using a WGS tumour-informed assay demonstrated high rates of MRD detection in patients with clinical stage I NSCLC
- In the absence of a direct comparison, this is higher than a previous report for a larger cohort using a first-generation WES-based assay<sup>1</sup>
- ctDNA-MRD detection was associated with RFS and provided clinically meaningful lead-times
- These results support exploring an ultrasensitive MRD-guided approach in stage I NSCLC to escalate adjuvant therapy



<sup>1</sup>Khan et al. ASCO 2024





# Thank You

- We thank all participating patients, their families, participating investigators and their site staff (University Health Network, Michael Garron, Unity Health)
- We gratefully acknowledge support from the Princess Margaret Cancer Foundation Grand Challenge, the Ontario Institute for Cancer Research CATA, and Foresight Diagnostics

